



Example 1 - LEG BAG STRAP

- Hook and loop sewn onto elastic base material




Example 2 - LEG BAG STRAP

- Hook-receptive elastic base material with silicone backing
- Sewn hook tip



Example 3 - CATHETER LEG BAG STRAP

- Hook-receptive elastic base material
- Welded hook tip



VELCRO® Brand

Urinary Bag Straps

A WIDE RANGE OF CUSTOMIZED OPTIONS
TO OFFER YOU THE BEST SOLUTION

Benefits

- Soft & comfortable materials ideal for skin contact
- Stretch fabrics that allow easy adjustment to the leg
- Variety of strengths to avoid circulatory problems and vein compression
- Silicone backings that avoid slipping of the bag
- Customization and branding

Options

- Colours
- Widths
- Different elasticity and elongations
- With backings (silicone)
- Cycle life
- Round or straight tips
- Easy peel feature
- Sewable / weldable strap assembly
- Customization with printing or sewn labels
- Packaging: in bulk, folded, packed in polybags



VELCRO® Brand

Catheter Securement **Concept**

PRACTICAL DEVICE DESIGNED TO ENSURE FIXATION OF CATHETER, TUBES OR LINES AND READJUSTMENT EVEN BY GLOVED PERSONNEL

Benefits

- Strong fixation of catheters, tubes and lines securely in place
- Soft & comfortable materials ideal for skin contact
- Simple application and handling even by gloved personnel
- VELCRO® Brand Hook & Loop closure allows readjustment
- Easy application onto patient skin by its adhesive backing
- Latex free
- Customization and branding

Options

- Shapes, dimensions and materials
- Different closure strengths
- Coloured or transparent closing tape
- Skin friendly adhesive alternatives: from 3 or 7 days; for infants and adults
- Adhesive resistant to incubators
- Customization with shapes or printing
- Packaging: individually, in bulk, in bags, own label


VELCRO® BRAND
HOOK & LOOP SECUREMENT

ANTI-SLIDDING ADHESIVE

SKIN FRIENDLY ADHESIVE



Example 1 - SHORT TERM ADHESIVE PATCH



Example 2 - LONG TERM ADHESIVE PATCH



Ask us about the VELCRO®
Brand Solutions for
direct skin contact as:

- Tube, cable and sensor holders



VELCRO® Brand

Surgical Tube **HOLDERS**

A WIDE RANGE OF CUSTOMIZED OPTIONS TO OFFER
YOU THE BEST SOLUTION

Benefits

- Designed for quick & secure placement of all sizes of tubing, lines and cords
- Fast removable of tubes by gloved personnel
- Re-adjustable hook & loop fastener
- Easy-to-apply adhesive base attaches firmly to skin, cloth or drapes
- Can be cut to the required length
- Customization and branding

Options


- Shapes and dimensions
- Strength materials and adhesives
- Skin friendly adhesive
- Laminate based to “soften” tubes
- Welded together or mated as sandwich
- Cycle life
- Round or straight tips
- Easy peel feature of adhesive (release liner) or strap (opening and closing)
- Sewable / weldable strap assembly
- Customization with printing or sewn labels
- Packaging: in bulk, folded, packed in polybags



Knit & extruded plastic hook tape

Example 1 - MATED SANDWICH CUT PIECE

- Knit and extruded plastic hook tape



Example 2 - WELDED STRAP

- Easy peel features



Ask us about the VELCRO®
Brand Solutions for...

- Surgical gowns closure

Example 3 - CUSTOMIZED SHAPE

- Soft laminate piece

